



For Anions: Use the KISS table. However, the second table has a test for nitrate.





Complexation



Selected Complexation Tests

p it simple science A "<u>co-ordination complex</u>" (or "metal complex") is a chemical compound or ion formed when a central metal atom becomes surrounded by, and bonded to, a number of molecules or ions. The surrounding molecules or ions are called "<u>ligands</u>" and may include polyatomic ions.

The "transition metals" (d-block of Periodic Table) often form complexes which have bright & characteristic colours. These complexes have long been used to make pigments for colouring fabrics, paints, etc.

They are also very useful for chemical identification...

Identifying Silver or Lead Ions

Silver ions form a white precipitate with Cl⁻ ions.

To verify that the precipitate IS AgCI (and not PbCI₂) carry out the precipitation in a test tube, then:

Add a few drops of dilute nitric acid to the pptte.
 Now add (drop-wise) dilute ammonia solution.
 (This is the same as ammonium hydroxide soln.)

Solid silver chloride will re-dissolve as it forms a soluble complex ion:

Lead chloride does NOT re-dissolve in ammonia.

To verify suspected Pb²⁺ ions use another similar complexation reaction:

1. Take a fresh sample of the "unknown". Add NaOH solution. This produces a pptte of Pb(OH) $_{2(s)}$.

$$Pb^{2+}_{(aq)} + 2OH^{-}_{(aq)} \longrightarrow Pb(OH)_{2(s)}$$

2. Continue adding NaOH solution. Lead(II) hydroxide will re-dissolve as it forms a soluble complex ion:

$$Pb(OH)_{2(s)} + 2OH^{-} \longrightarrow [Pb(OH)_{4}]^{2^{-}}_{(aq)}$$

Phosphate Ion Test

Phosphate ions can be positively identified as follows:

- In a test tube, add a few drops of nitric acid to your "unknown". Then add <u>ammonium molybdate</u> solution and heat gently in a bunsen flame.
- If phosphate ions are present a bright yellow precipitate will form.

Explanation:

Ammonium molybdate is a "metal complex" around an atom of the element molybdenum (atomic No. 42). Formula is $(NH_4)_2MoO_4$.

The yellow pptte is an even more complicated complex called ammonium phosphomolybdate. (The KISS Principle prevents us from going there!)

Identifying Iron Ions

Ferrous ions (Fe²⁺) and ferric ions (Fe³⁺) can be difficult to tell apart. This may be due to the rapid oxidation of Fe²⁺ to Fe³⁺ when solutions are exposed to air. This means that solutions containing ferrous ions soon contain ferric ions as well.

A definitive test for Fe³⁺ ions involves the formation of a complex ion with an intense red colour often described as "dragon's blood".

 Add a few drops of ammonium thiocyanate (NH₄CNS) to a fresh sample of your "unknown".

The bright red complex is soluble.

- The ferrous ion will NOT react. However, because it rapidly oxidises to ferric ion, you may get a palered colour in your test tube.
- There is another complexation test which is specific for ferrous ions. We leave it to your teacher, or your own research, to find it.

(The KISS Principle demands that we are NOT the ones to cause any confusion!)

Chloride, Bromide & Iodide Ions

All 3 of these will form a pptte with Ag⁺ ions:

Furthermore, if the precipitate is treated with ammonia (NH_3) the same complex ion is formed around the silver ion:

The difference is that the:

- chloride version is soluble, so the pptte re-dissolves.
- bromide version is "sparingly soluble" so it redissolves only if treated with concentrated NH₃.
- iodide version is highly insoluble. The pptte does not re-dissolve.

Now Try Worksheets 1 & 2